

SuperPette

Fully Autoclavable Pipettors

The new SCILOGEX SuperPette fully autoclavable digital pipettors have an excellent ergonomic and ultralight design providing effortless, ergonomic, accurate and precise pipetting experience. In addition, the volume lock helps protect from accidental volume change. The clear digital display is centrally located, for left and right handed users. UV and chemical resistant.



Smooth, Low Force, Effortless Button

Smooth Knob Adjustment

Light and smooth volume adjustment

Ergonomic Design

Comfortable & Effortless, Precise & Accurate

Pull up the lock wheel to set the volume

One Key Volume Lock Setting

Push down the lock wheel to lock the volume

Volume Display Window

New centrally located volume display

Simple Disassembly

Easy to maintain and calibrate



Ultra-light ergonomic design

Light weight and light pipetting force with soft range adjustment and piston movement helps prevent pipetting fatigue.



Accurate pipetting

Special material and the updated components & technologies effectively reduces the operating force for accuracy, precision and reliability.



Complete sterilization

Fully autoclavable and excellent protection against chemicals & UV light.

With Filters

Reduces volatile liquids entry into the pipette body; helps avoid liquid back-flow caused by improper operation and thus prevents piston corrosion and contamination



Volume Display Window

New centrally located volume display



Forward pipetting Suitable for low-viscosity liquids, press the plunger to the first stop when aspirating, and press the button to the second stop when dispensing.

Reverse pipetting Suitable for higher-viscosity liquids, easy to foam or micro samples, with higher pipetting accuracy. Press the plunger to the second stop when aspirating, and press the plunger to the first stop when dispensing.

Range selection Choose the appropriate range according to the pipetting volume. 30% - 100% of the maximum range of the pipette is the suitable range.

Specifications

Single-channel Adjustable Volume Pipettes (9 models, covering 0.1µl~10ml.)						
Volume Range	Increment	Test Volume	Accuracy Error		Precision Error	
			%	µL	%	µL
0.1~2.5µL	0.002µL	2.5µL	1.4%	0.04	0.7%	0.02
		1.25µL	2.5%	0.03	1.5%	0.02
		0.25µL	12.0%	0.03	6.0%	0.02
0.5~10µL	0.01µL	10µL	1.0%	0.10	0.4%	0.04
		5µL	1.5%	0.08	0.8%	0.04
		1µL	2.5%	0.03	1.8%	0.02
2~20µL	0.02µL	20µL	1.0%	0.20	0.3%	0.06
		10µL	1.2%	0.12	0.6%	0.06
		2µL	5.0%	0.10	1.5%	0.03
10~100µL	0.1µL	100µL	0.8%	0.80	0.2%	0.20
		50µL	1.0%	0.50	0.3%	0.15
		10µL	3.0%	0.30	1.0%	0.10
20~200µL	0.2µL	200µL	0.6%	1.20	0.2%	0.40
		100µL	1.0%	1.00	0.3%	0.30
		20µL	2.5%	0.50	0.7%	0.14
30~300µL	2µL	1000µL	0.6%	1.80	0.2%	0.60
		500µL	1.0%	1.50	0.3%	0.45
		100µL	2.5%	0.75	0.7%	0.21
100~1000µL	1µL	1000µL	0.6%	6.00	0.2%	2.00
		500µL	1%	5.00	0.2%	1.00
		100µL	3%	3.00	0.6%	0.60
1000~5000µL	5µL	5000µL	0.6%	30.00	0.2%	10.00
		2500µL	0.6%	15.00	0.3%	7.50
		1000µL	1.2%	12.00	0.3%	3.00
2000~10000µL	10µL	10000µL	0.6%	60.00	0.2%	20.00
		5000µL	0.8%	40.00	0.2%	10.00
		2000µL	1.5%	30.00	0.3%	6.00

* For accurate and reliable pipette calibration SuperPette conforms to ISO 8655:2002 standards.

