



## COVID-19 Response Solutions

### Rapidly Deployable Laboratory Space

In response to the unprecedented threat of the COVID-19 outbreak, Germfree offers the bioGO® line of turnkey Biocontainment Laboratories for immediate deployment.

- Applications include:
- Emergency Disease Surveillance
  - Added Capacity / Temporary Laboratory Space
  - Laboratory Swing Space
  - Any application where microbiological containment laboratory space is required



**53' BSL-2+ Mobile Biocontainment Laboratory**

- Drop Deck Dry Van Trailer Platform
- 5 Distinct Work Zones
  - BSL-2 Enhanced Laboratory
  - BSL-2 Laboratory (PCR Prep)
  - 2 Anterooms
  - Laboratory Command Center
- Suited for immediate deployment within the continental United States



**40' ISO Container BSL-3 Laboratory Building**

- High-Cube ISO Shipping Container Platform
- 3 Distinct Work Zones
  - WHO BSL-3 Laboratory
  - Anteroom
  - Main Entry Area
- Suited for immediate deployment worldwide



**24' BSL-2 Mobile Biocontainment Laboratory**

- Tag Trailer Platform
- 2 Distinct Work Zones
  - BSL-2 Laboratory
  - Anteroom
- Suited for immediate deployment within the continental United States



**BSL-3 Mobile Biocontainment Laboratory**

- 44' Trailer Platform
- 4 Distinct Work Zones
  - BSL-3 Laboratory
  - BSL-2 Laboratory
  - BSL-3 Anteroom/Shower Room
  - Laboratory Command Center
- Suited for immediate deployment worldwide

**+1 386.265.4300**  
**www.Germfree.com**