

Ultra Series

Vaccine & Pharmacy
Upright Glass Door
Refrigerators

Our **ULTRA CliniCool Series** medical refrigerators are specifically designed to provide the safest solution for storing vaccines and other medical cold storage products. These units provide exceptional temperature uniformity and security to ensure the efficacy & safety of your valuable samples. acceptable for use in non-

CliniCool
LABRepCo



LABRepCo

One Company • Many Solutions

101 Witmer Road, Suite 700 Horsham, PA 19044 / www.labrepco.com / (800) 521-0754 / sales@labrepco.com

Our **ULTRA Basic CliniCool Series** glass door refrigerators provide exceptional temperature uniformity for storing product in the medical applications market. These units come equipped with a microprocessor controller & digital temperature display, while featuring high/low temperature alarms. Available from 12 to 72 cu. ft. of internal storage capacity.

Standard Product Features:

- Cycle defrost
- Pharmacy/Vaccine Glass Door Refrigerator
- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- Adjustable operating temperature range: +36°F to +46°F (+2°C to +8°C)
- Alarm testing capabilities
- Freestanding thermometer with 3 year certificate of calibration
- 4 Adjustable shelves (adjustable in ½" increments)
- Magnetic door gasket for positive seal
- Heavy duty door handle
- 2" Swivel casters
- LED interior lights, safety shielded and switch controlled
- Probe access port (3/8")
- Powder coated steel white interior and exterior
- High density urethane foam cabinet insulation
- CFC-free refrigerant (R134a), foam insulation and packaging
- Pharmacy refrigerator/freezer toolkit & temperature logs
- Forced draft air circulation
- Warranty: 2 years parts and labor, plus an additional 3 years compressor parts only



Optional Features:

- Additional full shelf
- NIST Traceable® Thermometers
- NIST Traceable® Data Loggers

LHU-23-HG-PH

LABRepCo

One Company • Many Solutions

Specifications and Ordering Information:



Digital Temp Display/
Microprocessor Temp Controller



2" Swivel Casters



Models	LHU-12-HG-PH	LHU-16-HG-PH	LHU-23-HG-PH	LHU-26-HG-PH	LHU-33-SG-PH	LHU-36-HG-PH	LHU-47-SG-PH	LHU-49-HG-PH	LHU-69-SG-PH	LHU-72-HG-PH
Shipping Weight (lbs)	254	282	331	345	419	430	550	550	741	741
Overall Height (in.)	63 3/4	77	79 3/4	79 3/4	79 3/4	79 3/4	79 3/4	79 3/4	79 3/4	79 3/4
Overall Width (in.)	25	25	26 7/8	28 3/8	39 5/8	39 5/8	54	54	81	81
Overall Depth (in.)	29 3/4	29 3/4	34 3/4	36 1/2	33 1/4	34 3/4	33 1/4	34 3/4	33 1/4	34
Interior Cu. Ft.	12	16	23	26	33	36	47	49	69	72
Number of Shelves	3 adj	3 adj	4 adj	4 adj	8 adj	8 adj	8 adj	8 adj	12 adj	12
Defrost Type	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle
Refrigerant	R134a	R134a	R134a	R134a	R134a	R134a	R134a	R134a	R134a	R134a
Compressor Size HP	3/8	3/8	3/8	3/8	3/8	1/2	1/2	1/2	3/4	3/4
Plug (receptacle)	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R
Total Amp Draw	6	6	6	6	6	9	7	7	10	10
Voltage	115V, 3/8PH, 60Hz	115V, 3/8PH, 60Hz	115V, 1PH, 60Hz	115V, 3/8PH, 60Hz	115V, 3/8PH, 60Hz	115V, 3/8PH, 60Hz	115V, 1/2PH, 60Hz	115V, 1/2PH, 60Hz	115V, 3/4PH, 60Hz	115V, 3/4PH, 60Hz

Receptacle Required

115V/1PH/ 60Hz
NEMA-5-15R



Specifications subject to change without notice.



LABRepCo
One Company • Many Solutions

UltraElite

Vaccine & Pharmacy Upright Glass Door Refrigerators

Our **ULTRA ELITE CliniCool Series** glass door refrigerators provide the same exceptional temperature uniformity as our basic series, but include an upgraded digital temperature display and alarm module with additional features. These units are available from 12 to 72 cu. ft. of internal storage capacity.

Standard Product Features:

- Cycle defrost
- Microprocessor temperature controller
- Adjustable operating temperature range: +36°F to +46°F (+2°C to 8°C)
- Freestanding thermometer with 3 year certificate of calibration
- Standard Digital Temperature Display and Alarm Module
- Sample simulated digital temperature display; .1° Resolution
- Adjustable temperature display in Celsius or Fahrenheit
- Audible and visual high/low temperature alarms
- Sensor error alarm
- Power failure alarm
- Resettable 1-point min/max history
- Remote alarm contacts
- Alarm testing capabilities
- Alarm mute with ring back
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle
- 2" Swivel casters
- LED interior lights, safety shielded and switch controlled
- Probe access port (3/8")
- CFC-free refrigerant (R134a), foam insulation and packaging
- Pharmacy refrigerator/freezer toolkit & temperature logs
- Forced draft air circulation
- Warranty: 2 years parts and labor, plus an additional 3 years compressor parts only



LHE-23-HG-PH

Optional Features:

- Additional full shelf
- NIST Traceable® Thermometers
- NIST Traceable® Data Loggers

LABRepCo

One Company • Many Solutions

Specifications and Ordering Information:



Standard Digital Temp Display & Alarm Module




2" Swivel Casters



Models	LHE-12-HG-PH	LHE-16-HG-PH	LHE-23-HG-PH	LHE-26-HG-PH	LHE-33-SG-PH	LHE-47-SG-PH	LHE-49-HG-PH	LHE-69-SG-PH	LHE-72-HG-PH
Shipping Weight (lbs.)	258	275	336	355	419	550	538	741	743
Overall Height (in.)	63 3/4	77	79 3/4	79 3/4	79 3/4	79 3/4	79 3/4	79 3/4	79 3/4
Overall Width (in.)	25	25	26 7/8	28 3/8	39 5/8	54	54	81	81
Overall Depth (in.)	29 3/4	29 3/4	34 3/4	36 1/2	33 1/4	33 1/4	34 3/4	33 1/4	34
Interior Cu. Ft.	12	16	23	26	33	47	49	69	72
Number of Shelves	3 adj	3 adj	4 adj	4 adj	8 adj	8 adj	8 adj	12 adj	12 adj
Defrost Type	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle	Cycle
Refrigerant	R134a	R134a	R134a	R134a	R134a	R134a	R134a	R134a	R134a
Compressor Size HP	3/8	3/8	3/8	3/8	3/8	1/2	1/2	3/4	3/4
Plug (receptacle)	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R	5-15R
Total Amp Draw	6	6	6	6	6	7	7	10	10
Voltage	115V, 3/8PH, 60Hz	115V, 3/8PH, 60Hz	115V, 1PH, 60Hz	115V, 3/8PH, 60Hz	115V, 3/8PH, 60Hz	115V, 3/8PH, 60Hz	115V, 1/2PH, 60Hz	115V, 3/4PH, 60Hz	115V, 3/4PH, 60Hz

Specifications subject to change without notice.

Receptacle Required
115V/1PH/ 60Hz
NEMA-5-15R




LABRepCo
One Company • Many Solutions