



**Nexxit Chair 9033-A,  
Desk Height,  
Freewheeling Casters,  
Black Vinyl, With 4D  
Arms**

**Article #: LRC-NXT33A\***

*\*In-Stock Quick Ship Item*

---

**Features**

The Nexxit is optimized for use in lab environments with a modern design, ergonomic comfort for long sit times and sculpted shoulder areas to assist with a researchers' range of motion. The light scale offers ease of movement and full adjustability to each user while engineered to stand up to the numerous rigors within the lab.

---

**Specifications**

<b>Height Adjustment:</b>	Standard Desk Height Range 17-3/4" - 23-5/8"
<b>Seat Depth Range:</b>	16-5/8" - 18-7/8"
<b>Handle Color:</b>	In-Stock (Blue); Factory Orders (Orange, Green and Grey)
<b>Frame Color:</b>	Black
<b>Seat Material:</b>	Black Vinyl
<b>Castors:</b>	Freewheeling
<b>Base:</b>	Black Plastic Base
<b>Cleanroom Rating:</b>	ISO Level 7
<b>Bio-Safety BSL:</b>	Level 2
<b>BIFMA Rating:</b>	Complies with Level e3 Platinum
<b>SEFA Rating:</b>	Complies with Scientific Equipment & Furniture Assn. - SEFA12
<b>Weight Capacity:</b>	300 lbs.
<b>Product Weight:</b>	33 lbs.
<b>Warranty:</b>	10-Year Warranty
<b>Original Design:</b>	German Engineering and Production
<b>ESD Capable:</b>	Available by Request