



**Nexxit Chair Plus 9033+,
Desk-to-Bench Height,
Freewheeling Casters,
Black Vinyl, Foot
Ring, With 4D Arms**

Article #: LRC-NXT33PA*

**In-Stock Quick Ship Item*

Features

The Nexxit is optimized for use in lab environments with a modern design, ergonomic comfort for long sit times and sculpted shoulder areas to assist with a researchers' range of motion. The light scale offers ease of movement and full adjustability to each user while engineered to stand up to the numerous rigors within the lab.

Specifications

Seat Adjustment:	Extended Desk-to-Bench Height Range 20-1/2" - 28"
Seat Depth Range:	16-5/8" - 18-7/8"
Handle Color:	In-Stock (Blue); Factory Orders (Orange, Green and Grey)
Frame Color:	Black
Seat Material:	Black Vinyl
Castors:	Freewheeling
Base:	Black Plastic Base
Cleanroom Rating:	ISO Level 7
Bio-Safety BSL:	Level 2
BIFMA Rating:	Complies with Level e3 Platinum
SEFA Rating:	Complies with Scientific Equipment & Furniture Assn.- SEFA12
Weight Capacity:	300 lbs.
Product Weight:	33 lbs.
Warranty:	10-Year Warranty
Original Design:	German Engineering and Production
ESD Capable:	Available by Request