



Lab Vortex Mixers Accessories

Catalog



***Explore the various accessories
that you can use with your vortex mixer
Holders, Microplates, Foam Inserts, and more!***

Ingeniously Practical

WWW.OHAUS.COM

VORTEX MIXER ACCESSORIES

Micro-Tube Holder

Mixes (48) 0.25 to 2 mL micro-tubes.
Requires Insert Retainer.

Description	Item No.
Micro-Tube Holder (2 pack)	30400232
Insert Retainer	30400227

9 to 13 mm Tube Holder

Ideal for mixing 5 mL culture tubes and micro-vials.
Requires Insert Retainer.

Description	Item No.
9 to 13 mm Tube Holder (2 pack)	30400229
Insert Retainer	30400227

14 to 19 mm Tube Holder

Ideal for mixing up to (8) 15 mL centrifuge tubes.
Requires Insert Retainer.

Description	Item No.
14 to 19 mm Tube Holder (2 pack)	30400230
Insert Retainer	30400227

20 to 25 mm Tube Holder

Ideal for mixing up to (8) 50 mL centrifuge tubes.
Requires Insert Retainer

Description	Item No.
20 to 25 mm Tube Holder (2 pack)	30400231
Insert Retainer	30400227

Vessel Harness

Mixes Erlenmeyer flasks and media bottles.
Requires Insert Retainer.

Description	Item No.
Vessel Harness (2 pack)	30400228
Insert Retainer	30400227

Cup Head

Designed for mixing 1 tube at a time.

Description	Item No.
Cup Head	30400235

3" Rubber Head Cover and 3" Head

Designed for mixing irregular shaped objects.

Description	Item No.
3" (7.6 cm) Rubber Head Cover	30400196
3" (7.6 cm) Head	30400236
3" (7.6 cm) Head with Rubber Cover	30400237

Microplate Holder

Ideal for mixing 96-well plates or deep well blocks.
Requires Insert Retainer.

Description	Item No.
Microplate Holder (2 pack)	30400233
Insert Retainer	30400227

Flat Foam Insert

Ideal for custom applications. Can be cut or drilled to fit your specifications. Requires Insert Retainer.

Description	Item No.
Flat Foam Insert (2 pack)	30400234
Insert Retainer	30400227

Microwell Kit with Retainer

Includes:
• 2 Microplate Holders
• 1 Insert Retainer

Description	Item No.
Microwell Kit with Retainer	30400199

Microtube Kit with /Retainer

Includes:
• 2 Micro-Tube Holders
• 1 Insert Retainer

Description	Item No.
Microtube Kit with /Retainer	30400200

Foam Insert Variety Pack

Includes:
• (1) 9 to 13 mm Tube Holder
• (1) 14 to 19 mm Tube Holder
• (1) 20 to 25 mm Tube Holder
• 1 Flat Foam Insert

Description	Item No.
Foam Insert Variety Pack	30400201

Micro-Tube and Microplate Holder Kit

Includes:
• 1 Micro-Tube Holder
• 1 Microplate Holder
• 1 Insert Retainer

Description	Item No.
Micro-Tube & Microplate Holder Kit	30400198

Tube Holder Kit

Includes:
• (1) 9 to 13 mm Tube Holder
• (1) 14 to 19 mm Tube Holder
• (1) 20 to 25 mm Tube Holder
• 1 Flat Foam Insert
• 2 Vessel Harnesses
• 1 Insert Retainer

Description	Item No.
Tube Holder Kit	30400197

OHAUS Vortex Mixer accessories are interchangeable with other manufacturers' models. Call 1-800-672-7722 for assistance.

VORTEX MIXER ACCESSORIES

Single Tube Holder

Single tube, hands free mixing. Easily attaches to the top of any Vortex Mixer with the use of a strong magnetic base. Accepts tubes from 2.5 to 4.5" (6.4 to 11.4 cm) in length. Minimum tube diameter of 0.75" (19 mm).

Description	Item No.
Single Tube Holder	30400202

Single Tube Holder Accessories

(requires Single Tube Holder)

0.5 mL Micro-Tube Holder

Mixes (24) 0.5 mL micro-tubes.
For use with Single Tube Holder.

Description	Item No.
Micro-Tube Holder (0.5 mL)	30400204

Ampule Tube Holder

Mixes up to 4 storage vials and test tubes.

Description	Item No.
15 to 17 mm Ampule Tube Holder	30400207
10 to 17 mm Ampule Tube Holder	30400208

Stainless Steel Microtube Holder

Mixes up to (12) 1.5 to 5 mL tubes.
Stainless steel construction.

Description	Item No.
Microtube Holder	30400206

Adapter for Vortex-Genie® Mixer*

Adapter plate easily adheres to the Vortex-Genie® Mixer housing so Single Tube Holder (sold separately) can be attached.

Description	Item No.
Adapter plate	30400203

*The Vortex-Genie® Mixer is a registered trademark of Scientific Industries, Inc.

1.5 mL to 2.0 mL Micro-Tube Holder

Mixes (18) 1.5 to 2 mL micro-tubes.
For use with Single Tube Holder.

Description	Item No.
Microtube Holder (1.5 mL to 2.0 mL)	30400205

HEAVY DUTY AND MICROPLATE VORTEX MIXER ACCESSORIES

Foam Insert for 0.5 mL Microtubes

Foam insert holds (52) 0.5 mL microtubes.
Requires Universal Holder.

Description	Item No.
0.5 mL microtubes	30400216
Universal Holder	30400226

Foam Insert for 1.5 to 2.0 mL Microtubes

Foam insert holds (38) 1.5 to 2.0 mL microtubes.
Requires Universal Holder.

Description	Item No.
1.5 to 2.0 mL microtubes	30400217
Universal Holder	30400226

Foam Insert for 12-13 mm Test Tubes

Foam insert holds (34) 12-13 mm diameter test tubes.
Requires Universal Holder.

Description	Item No.
12-13 mm test tubes	30400220
Universal Holder	30400226

Foam Insert for 15-18 mm Test Tubes

Foam insert holds (20) 15-18 mm diameter test tubes.
Ideal for 15 mL centrifuge tubes.
Requires Universal Holder.

Description	Item No.
15-18 mm test tubes	30400221
Universal Holder	30400226

Foam Insert for 19-21 mm Test Tubes

Foam insert holds (18) 19-21 mm diameter test tubes.
Requires Universal Holder.

Description	Item No.
19-21 mm test tubes	30400222
Universal Holder	30400226

Foam Insert for 22-25 mm Test Tubes

Foam insert holds (13) 22-25 mm diameter test tubes.
Requires Universal Holder.

Description	Item No.
22-25 mm test tubes	30400223
Universal Holder	30400226

Foam Insert for 26-29 mm Test Tubes

Foam insert holds (4) 26-29 mm diameter test tubes. Ideal for 50 mL centrifuge tubes. Requires Universal Holder.

Description	Item No.
26-29 mm test tubes	30400224
Universal Holder	30400226

Cup Head*

Designed for mixing 1 tube at a time.

Description	Item No.
Cup Head	30400210

Small Vessel Holder

Rubber holder secures 125 and 250 mL Erlenmeyer flasks.
Vessel holder also includes a grip mat.
Requires Universal Holder.

Description	Item No.
Small Vessel Holder	30400218
Universal Holder	30400226

Large Vessel Holder

Rubber holder secures 500 and 1000 mL Erlenmeyer flasks.
Vessel holder also includes a grip mat.
Requires Universal Holder.

Description	Item No.
Large Vessel Holder	30400211
Universal Holder	30400226

Microplate Holder (Single)*

Designed to hold one standard microplate.

Description	Item No.
Single Holder	30400215

Microplate Holder (Double)*

Designed to hold two standard microplates.

Description	Item No.
Double Holder	30400213

Microplate Holder (Quad)

Designed to hold four standard microplates.

Description	Item No.
Quad Holder	30400214

Stackable Microplate Holder Four

Designed to maximize the capacity of the Heavy-Duty Vortex Mixer to eight microplates by stacking the tray on top of the four plate holder.

Description	Item No.
Stackable Holder	30400212

Flat Foam Insert

Ideal for custom applications. Can be cut or drilled to fit your specifications. Requires Universal Holder.

Description	Item No.
Flat Foam Insert	30400209
Universal Holder	30400226

Universal Holder & Cover


Replacement for items supplied with Heavy-Duty Vortex Mixer. Cover allows for mixing irregularly shaped objects.

Description	Item No.
Universal Holder	30400226
Universal Holder Cover	30400225

Unless noted with an asterisk (), for use on Heavy-Duty Vortex only.

MULTI-TUBE VORTEXER ACCESSORIES

Foam Test Tube Racks



Description	Tube Capacity	Color	Dimensions (L x W x H)	Item No.
10 mm Test Tube Foam Rack	50	Gray	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400238
12 mm Test Tube Foam Rack	50	Blue	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400239
13 mm Test Tube Foam Rack	50	Yellow	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400240
16 mm Test Tube Foam Rack (for 15 mL centrifuge tubes)	50	Green	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400241
25 mm Test Tube Foam Rack	28	White	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400243
29 mm Test Tube Foam Rack (for 50 mL centrifuge tubes)	15	Red	5.5 x 9.5 x 2" (14 x 24.1 x 5.1 cm)	30400242
Replacement Tray Pad Set (upper & lower)	N/A	Gray	7 x 12 x 1" (17.8 x 30.5 x 2.5 cm)	30400245

Post Extension Kit

Adds 6" (15.2 cm) to the post of the Multi-Tube Vortex Mixer to accommodate tubes up to 10" (25.4 cm) tall.


Description	Item No.
Post Extension Kit	30400244





Visit the Lab Vortex Accessories Product Page:

Click the image or use the QR code



Contact Us:

OHAUS CORPORATION

Tel: 800.672.7722
 973.377.9000
 Fax: 973.944.7177

OHAUS Europe GmbH

e-mail: ssc@ohaus.com
 Tel: 0041 22 567 53 19
 e-mail: tsc@ohaus.com
 Tel: 0041 22 567 53 20

OHAUS Asia/Pacific Headquarters

e-mail: APmarketing@ohaus.com
 Tel: +86 21 64855408

奥豪斯仪器 OHAUS Instruments

客户咨询邮箱: inquirychina@ohaus.com
 售后服务邮箱: ohauservice@ohaus.com

www.ohaus.com

80776906_A 20240229 © Copyright OHAUS Corporation

The management system governing the manufacture of this product is ISO 9001:2015 certified.

